

Portable pH Meter – MSpH 816



MS pH 816 is a compact Handy highly precision pH meter for continuous 0-14 pH reading & is easy to operate as it is having least number of controls. Auto Bat. Low Indication is displayed as BAT when battery voltage falls to 7.6 volts dc.

Specifications :-

pH Mode

Range	0.00 to 14.00 pH
Resolution	0.01 pH
Accuracy	0.05 pH \pm 1 digit
Input Impedance	5×10^{11} Ohms
Temperature Control	0 - 80°C Manual
Display	3.5 digit LCD display
Calibration	Two buffers calibration (manually) 7pH & 4pH Or 9.2pH
Power Requirement	9V dry cell flat pack Battery
Environment	5 to 50c
Dimension	78x105x40mm (ABS Cabinet)
Weight	500gms (Approx)

Note : Due to Continues improvement in Product the Specifications may vary