


Table Top pH Meter MS pH 616

*Table Top pH meter
with pH electrode and stand*



It is a high-precision pH meter for continuous 0-14 pH reading & milli volts. The pH meter is easy to operate as it has the least number of controls. CHECK facility saves time from every day buffer standardization and hence is most suitable for quick analysis. It is a very sleek attractive finished instrument.

Specifications:-

	pH Mode	mV Mode
Range	0.00 to 14.00 pH	0 to ± 1999 mV
Resolution	0.01 pH	± 1 mV
Accuracy	0.01 pH ± 1 digit	1 mV ± 1 digit
Input Impedance	10^{12} Ohms	
Temperature Control	0 to 100°C Manual	
Display	3.5 digit LED display with auto polarity & decimal	
Calibration	Two buffers calibration (manually) 7pH & 4pH Or 9.2pH	
Power Requirement	230V A.C $\pm 10\%$, 50Hz single phase	
Environment	5 to 50°C	
Dimensions	205 x 65 x 130mm (Aluminum powder coated cabinet)	
Weight	1.5kgs (Approx) including accessories	
Accessories	pH Sensor, pH Sensor Stand, Dust Cover, Instruction Manual	